

## EXAMPLE A

# Use of conjugated Fab fragments for labeling and blocking

Read more about  
Fab fragments at:  
[jacksonimmuno.com](http://jacksonimmuno.com)

## Detection of two unlabeled primary antibodies from the same host species

Example A shows one of the possible protocols used for double labeling two unconjugated primary antibodies from the same host species. For more protocols visit: [www.jacksonimmuno.com/technical/products/protocols/double-labeling-same-species-primary](http://www.jacksonimmuno.com/technical/products/protocols/double-labeling-same-species-primary).

The success of these experimental designs will require some empirical manipulations. Optimizing reagent concentrations in each step or switching the labeling sequence of the two antigens may influence the outcome.

- Labeling the less abundant primary antibody first increases blocking efficiency.
- Blocking with an appropriate normal serum helps to reduce background.
- To avoid displacement of the Fab antibody by the labeled secondary antibody, a light post-fixation with glutaraldehyde may be used, provided that it does not affect antigenicity of the target proteins.


**Important note:** The monovalent Fab fragments have not been adsorbed to remove cross-reactivities to other species. If the experimental sample contains endogenous immunoglobulins [Example C](#) should be used. Example A or B could introduce background.




Jackson ImmunoResearch Laboratories Inc. is certified by BSI to ISO 9001:2015 under certificate number FM 545248.

## EXAMPLE A


Detection of two unlabeled primary antibodies from the same host species




1. After blocking with normal serum, incubate with the first primary antibody, in this example Rabbit Anti-Antigen X. Wash.



2. Incubate with excess conjugated secondary antibody, in this example Alexa Fluor® 488-Fab fragment Goat Anti-Rabbit IgG (H+L). Wash.



3. Incubate with the second primary antibody, Rabbit Anti-Antigen Y.



4. Incubate with a second conjugated secondary antibody, in this example Rhodamine Red™-X-Goat Anti-Rabbit IgG (H+L). Wash.

**Application notes:** (1) Monovalent Fab fragments have not been adsorbed against other species, so they may cross-react with endogenous Ig. Use Example C to avoid detection of endogenous Ig.

(2) Example A may require a high concentration of conjugated Fab to saturate the first primary antibody. If this results in unacceptable background, try a lower concentration of the conjugated Fab, followed by further blocking with unconjugated Fab.

## Key



Rabbit  
Anti-  
Antigen X



Rabbit  
Anti-  
Antigen Y



Goat  
Anti-Rabbit  
IgG (H+L)



Fab fragment  
Goat Anti-  
Rabbit IgG  
(H+L)



Alexa  
Fluor.  
488



Rhodamine  
Red-X



Antigen X



Antigen Y